




ATRAUMATIC  
EXTRACTION







## ANAQA ELEVATORS

- Unique 14 mm handle for best grip and handling.
- Circular handle improves handling at any position.
- Innovative etched design provides excellent grip.
- The depth and angles of the handle's grooves are optimised to provide the maximum friction between the instrument and the glove, virtually eliminating the risk of slippage during dental procedures.
- Made from high quality medical grade stainless steel with optional titanium carbide coated satin black finish, or luxurious titanium nitride coated gold finish.
- Highly efficient, sharp and slim tips for better access and function.

## ANAQA PERIOTOMES

- Wide range of Periotomes.
- Unique Periotomes with Flexible Tips made from high quality spring stainless steel with specialized heat treatment giving strength and flexibility.
- Interchangeable Periotomes - An economical choice.
- Standard Periotome set for more convenient use.
- Improved handling because of unique design of knurls and slightly increased length of the grip.
- The depth and angles of the handle's grooves are optimised to provide the maximum friction between the instrument and the glove, virtually eliminating the risk of slippage during dental procedures.
- Three different finishes because we care about aesthetics too.



## QUALITY YOU CAN COUNT ON

At ISFY, we don't take any chances when it comes to quality. That's why every Anaqa product is subjected to our rigorous quality control process.

We are so confident in our instruments that we offer an industry leading warranty of between ONE and FIVE years on all Anaqa products, for your total peace of mind.



**ONE YEAR WARRANTY**

for details please read the warranty policy



Luxation Instrument 3C  
(3mm curved)

210250SS	Standard Stainless Steel
210250NB	Nano Coating Black
210250NG	Nano Coating Gold



Luxation Instrument 3ST  
(3mm straight)

210251SS	Standard Stainless Steel
210251NB	Nano Coating Black
210251NG	Nano Coating Gold



Luxation Instrument 5C  
(5mm curved)

210252SS	Standard Stainless Steel
210252NB	Nano Coating Black
210252NG	Nano Coating Gold



Luxation Instrument 5ST  
(5mm straight)

210253SS	Standard Stainless Steel
210253NB	Nano Coating Black
210253NG	Nano Coating Gold



Luxation Instrument 2C  
(2mm curved)

210262SS	Standard Stainless Steel
210262NB	Nano Coating Black
210262NG	Nano Coating Gold



Luxation Instrument 2ST  
(2mm straight)

210263SS	Standard Stainless Steel
210263NB	Nano Coating Black
210263NG	Nano Coating Gold



Luxation Instrument 4C  
(4mm curved)

210264SS	Standard Stainless Steel
210264NB	Nano Coating Black
210264NG	Nano Coating Gold



Luxation Instrument 4ST  
(4mm straight)

210265SS	Standard Stainless Steel
210265NB	Nano Coating Black
210265NG	Nano Coating Gold







Warwick James Elevator  
(straight)

210254SS	Standard Stainless Steel
210254NB	Nano Coating Black
210254NG	Nano Coating Gold



Warwick James Elevator  
(left)

210255SS	Standard Stainless Steel
210255NB	Nano Coating Black
210255NG	Nano Coating Gold



Warwick James Elevator  
(right)

210256SS	Standard Stainless Steel
210256NB	Nano Coating Black
210256NG	Nano Coating Gold



Couplands Bone Chisel No. 1

210257SS	Standard Stainless Steel
210257NB	Nano Coating Black
210257NG	Nano Coating Gold



Couplands Bone Chisel No. 2

210258SS	Standard Stainless Steel
210258NB	Nano Coating Black
210258NG	Nano Coating Gold



Couplands Bone Chisel No. 3

210259SS	Standard Stainless Steel
210259NB	Nano Coating Black
210259NG	Nano Coating Gold



Cryers Elevator No. 30  
(left)

210260SS	Standard Stainless Steel
210260NB	Nano Coating Black
210260NG	Nano Coating Gold



Cryers Elevator No. 31  
(right)

210261SS	Standard Stainless Steel
210261NB	Nano Coating Black
210261NG	Nano Coating Gold







Root Tip Elevator  
(right)

210270SS	Standard Stainless Steel
210270NB	Nano Coating Black
210270NG	Nano Coating Gold



Root Tip Elevator  
(left)

210271SS	Standard Stainless Steel
210271NB	Nano Coating Black
210271NG	Nano Coating Gold



Root Tip Elevator  
(straight)

210272SS	Standard Stainless Steel
210272NB	Nano Coating Black
210272NG	Nano Coating Gold



Heidebrink Root Pick Elevator  
E13/E14

210281SS	Standard Stainless Steel
210281NB	Nano Coating Black
210281NG	Nano Coating Gold



Millers Bone File No. 52

210282SS	Standard Stainless Steel
210282NB	Nano Coating Black
210282NG	Nano Coating Gold



Periosteal Elevator Prichard 3

210278SS	Standard Stainless Steel
210278NB	Nano Coating Black
210278NG	Nano Coating Gold



Periosteal Elevator Buser

210279SS	Standard Stainless Steel
210279NB	Nano Coating Black
210279NG	Nano Coating Gold



Periodontal Chisel Rhodes  
36/37

210280SS	Standard Stainless Steel
210280NB	Nano Coating Black
210280NG	Nano Coating Gold





Periosteal Elevator Molt 9

210230SS	Standard Stainless Steel
210230NB	Nano Coating Black
210230NG	Nano Coating Gold



Pertotome 1

210283SS	Standard Stainless Steel
210283NB	Nano Coating Black
210283NG	Nano Coating Gold



Periotome 2

210284SS	Standard Stainless Steel
210284NB	Nano Coating Black
210284NG	Nano Coating Gold



Periotome 3

210285SS	Standard Stainless Steel
210285NB	Nano Coating Black
210285NG	Nano Coating Gold



Interchangeable Periotome Handle

210294SS	Standard Stainless Steel
210294NB	Nano Coating Black
210294NG	Nano Coating Gold

210320 - Interchangeable Periotome Set



Interchangeable Periotome Blade Wide

210295SS	Standard Stainless Steel
210295NB	Nano Coating Black
210295NG	Nano Coating Gold



Interchangeable Periotome Blade Narrow

210296SS	Standard Stainless Steel
210296NB	Nano Coating Black
210296NG	Nano Coating Gold



Interchangeable Periotome Blade Angled

210297SS	Standard Stainless Steel
210297NB	Nano Coating Black
210297NG	Nano Coating Gold







Flexible Periotome 2.5 mm  
(Right)

210286SS	Standard Stainless Steel
210286NB	Nano Coating Black
210286NG	Nano Coating Gold



Flexible Periotome 2.5 mm  
(Left)

210287SS	Standard Stainless Steel
210287NB	Nano Coating Black
210287NG	Nano Coating Gold



Flexible Periotome 2 mm  
(Curved)

210288SS	Standard Stainless Steel
210288NB	Nano Coating Black
210288NG	Nano Coating Gold



Flexible Periotome 3 mm  
(Curved)

210289SS	Standard Stainless Steel
210289NB	Nano Coating Black
210289NG	Nano Coating Gold



Flexible Periotome 2 mm  
(Straight)

210290SS	Standard Stainless Steel
210290NB	Nano Coating Black
210290NG	Nano Coating Gold



Flexible Periotome 3 mm  
(Straight)

210291SS	Standard Stainless Steel
210291NB	Nano Coating Black
210291NG	Nano Coating Gold



Flexible Periotome 2-3 mm  
(Lateral)

210292SS	Standard Stainless Steel
210292NB	Nano Coating Black
210292NG	Nano Coating Gold



Flexible Periotome 2-3 mm  
(Anterior)

210293SS	Standard Stainless Steel
210293NB	Nano Coating Black
210293NG	Nano Coating Gold



### Dental Syringe 2.2 ml (Self Aspirating - Slitting Type)

220351SS	Standard Stainless Steel
220351NB	Nano Coating Black
<b>B</b>	<b>RB</b>
<b>G</b>	<b>P</b>
<b>PN</b>	<b>RU</b>

Prefix Catalogue No. with the Colour Code



### Dental Syringe 1.8 ml (Self Aspirating - Slitting Type)

220352SS	Standard Stainless Steel
220352NB	Nano Coating Black
<b>B</b>	<b>RB</b>
<b>G</b>	<b>P</b>
<b>PN</b>	<b>RU</b>

Prefix Catalogue No. with the Colour Code

We develop our instruments through a continuous process of monitoring and responding to feedback from our customers, therefore our products are continually evolving. Our range of instruments cover most clinical situations, and all of our products are developed with performance and value at the core of our thought process. To discover the full range visit [www.isfy.co.uk](http://www.isfy.co.uk).







**Cleaning:** Clean and remove all debris immediately after use. Use an enzymatic cleaning solution in an ultrasonic cleaner to clean instruments. Your ultrasonic solution should be changed daily. Hinged instruments should be cleaned and sterilized in the open position. Discoloration and oxidation may occur if instruments are improperly cleaned. To increase the life of your instruments, clean them well. Use detergents and disinfectants suitable for use with medical and/or dental instruments.

**Lubrication:** To increase the life of your instruments, routine lubrication is recommended. Prior to sterilisation, lubricate instruments. This should be done on daily basis or at least weekly. Use only lubricants designed for precision hinged instruments. If using a dry heat steriliser, be sure to use a lubricant that is compatible with the operating temperatures of your dry heat unit.

**Sharpening:** Regular sharpening to your cutters will increase its service life, and its ability to cut more effectively. Frequency of sharpening depends on the frequency of use, as well as, the types and diameters of wires on which the cutter is being used. Check cutters frequently to assess if sharpening is needed.

**Sterilising:** Instruments should be in the open position when sterilised. The most efficient method is steam autoclave. Anaqa Range of instruments are actually designed for use with steam autoclave as well as Ethylene Oxide (EO) sterilisation. Steam autoclaves will not dull cutters. Other types of sterilising methods include chemical vapour or dry heat. Using cold sterilants is not recommended. If used, they can chemically attack your instruments.

**Corrosion:** Anaqa range of instruments are very resistant to corrosion. You can prevent corrosion problems by properly following all manufacturer's recommendations. Discoloration may appear that may not be corrosion. Organic materials that remain can give the appearance of corrosion and be mistaken for rust. An Orange/Brown stain is usually a phosphate layer caused by water sources, cleaning detergents, sterilisation solutions. Black stain is usually acid reaction caused by detergents. Dark brown stain usually caused by dried blood residues. Bluish / black stain usually is plating caused by autoclave cleaning material. Never mix a corroded instrument with new instruments as oxidation may be spread to new instruments. It is recommended that corroded instrument be discarded or replaced.

**Maintenance and Care:** It is extremely important to follow proper handling, care and sterilisation methods to ensure proper instrument function and improve the instrument life. Avoid tap water during the cleaning and sterilisation of instruments. Tap water contains high levels of minerals that can cause instrument discolouration and corrosion. Tap water can also neutralize rust inhibiting and lubricating solution. Always use purified, distilled water in ultrasonic cleaning and autoclave sterilisation units. It is also recommended to use a no-rinse ultrasonic solution that contains both a rust inhibitor and a lubricant prior to sterilisation and follow the manufacturer's exact solution measurements.

Do not leave instruments in the ultrasonic solution after cycle is complete, remove immediately and dry completely. Autoclaves contain high levels of moisture and can be damaging if not functioning properly. Instrument handling is also important. Make sure instruments are dry before sterilisation, keep jaws open, and do not overload trays. When cycle is complete, remove instruments immediately after cooling. Always check that joints and tips are moisture free or corrosion can occur. Lubricate with silicone-based lubricant before storage.

**Corrosion Prevention:** Corrosion can be an issue with any stainless steel instrument. Stainless steel requires Oxygen to form a protective chromium oxide surface layer, which protects against corrosion. Water droplets, organic residue and bonding agents left on the instruments prevent oxygen contact with the surface of instrument and protective layer of chromium oxide will not form. It is recommended

- Instruments should be thoroughly cleaned of all residual matter prior to sterilisation.
- All instruments should be sterilized in their open position.
- All instruments should be thoroughly air or towel dried taking special care in the joint areas and crevices.
- Instruments that have corrosion should never be sterilised with non-corroded instruments as the iron oxide on the corroded instrument can transfer to the non-corroded instruments permanently.



5 YEARS	Extraction Forceps
5 YEARS	Orthodontic Pliers & Cutters
1 YEAR	All Periodontic, Hygiene and Plastic Instruments

ISFY guarantees against defects in materials and manufacturing. If your Anaqa instruments should fail or are damaged due to a manufacturing or workmanship defect, we will repair or replace the product without quibble or charge once the faulty product is returned with proof of purchase.

Please note that this guarantee does not cover damage caused by accident, improper care, misuse, alterations to the Anaqa instruments or user negligence.

Anaqa extraction instruments should be disinfected following the processing guidelines provided. This guarantee does not cover damages which arise because these guidelines have not been followed.

ISFY reserves the right to decide whether an instrument is to be repaired or replaced. Where instruments require sharpening, this is considered routine maintenance with normal use and is not covered by any guarantee.

## WARRANTY INFORMATION

For more details, please visit [www.isfy.co.uk](http://www.isfy.co.uk) or contact us via email on [contact@isfy.co.uk](mailto:contact@isfy.co.uk).





# ISFY

QUALITY THROUGH INNOVATION

ISFY manufactures high quality, durable single and reusable surgical instruments. We also manufacture equipment to clients' exacting specifications.

We are serious about standards, and for that reason we have sourced the very best materials, machinery and processes from across the globe. For example, we use steel from Japan and tungsten carbide tips from the UK, for our premium Anaqa range.

We also invest in the very best people – from the highly skilled and experienced teams that manufacture our equipment, to the management team of the company, who between them have a wealth of experience within the industry.

We will never compromise on quality, which is why we offer industry-leading warranties on our products – up to FIVE years on our Anaqa range. Service is at the heart of everything we do. We pride ourselves on our flexibility and personal touch.

We are constantly looking to improve our processes and explore the latest technologies, materials and techniques to generate what we call “quality through innovation”.

ISFY Ltd  
I-Centre, Oakham Business Park  
Hamilton Way  
Mansfield  
Nottinghamshire  
NG18 5BR  
UK  
Tel: +44 (0)115 786 7111  
Email: [contact@isfy.co.uk](mailto:contact@isfy.co.uk)  
Web: [www.isfy.co.uk](http://www.isfy.co.uk)